

Nyhet! WelchAllyn® Blodtrycksmätare

767-Series Wall and Mobile Aneroid



- Certified accurate to ± 3 mm Hg.
- Lifetime calibration warranty.
- Jeweled movement for long life.
- Laser-engraved dial face for accuracy.
- Gear-free aneroid movement for reliability.
- Heat treated, nickel-beryllium diaphragm for lifetime stability and repeatability.
- Recessed dial face reduces parallax.
- Unique FlexiPort™ blood pressure cuff.
- Mercury-free and latex-free for safety.
- Large, easy-to-read dial face (15 cm).
- Built-in cuff storage basket on wall model.
- Mobile stand swivels 360° for clear viewing.



WelchAllyn®
Advancing Frontline Care™

Welch Allyn 767-Serien finns som vägg-, bords-, golvstativ-, och rälsmodell med FlexiPort™-, AJM™- eller TriCUFF®-manschett



FlexiPort™ familjen